



AGH UNIVERSITY OF SCIENCE
AND TECHNOLOGY

Optimisation of the X-ray fluorescence imaging system for mapping of pigments in historical paintings


Alicja Zielińska, B. Mindur, W. Dąbrowski

Faculty of Physics and Applied Computer Science
AGH University of Science and Technology


Jagiellonian Symposium of Fundamental and Applied Subatomic Physics
Kraków, 7 - 13 June 2015

Work supported by
The National Centre for Research and Development grant no. PBS3/A9/29/2015.

- 1 Introduction
 - ▶ Brief description of the motivation of our work
- 2 Detection system
 - ▶ Short overview of the Gas Electron Multiplier (GEM) based imaging system
- 3 Pigment specific imaging results
 - ▶ Presentation of the results obtained during scanning of painting phantom
- 4 Position resolution of the system
 - ▶ Information about the position resolution of the system in a few configurations
- 5 Conclusions
 - ▶ Remarks on results and a few words about future plans




Imaging of $10 \times 10 \text{ cm}^2$ area at fixed position!



The main advantages of the presented system are as follows

- Imaging of large areas of paintings in a fixed detector position (safety of the investigated objects)
- Relatively fast measurement (hours)
- Infinite depth of field (of camera obscura)
- Affordable costs

The principle of the method has already been proven.




Imaging of $10 \times 10 \text{ cm}^2$
area at fixed position!

Key components of the system

- Water-cooled X-ray tube with molybdenum or copper anode working at 20-30 kV
- Beam cross-section area $4 \times 10 \text{ cm}^2$ which allow to illuminate the painting area of about $10 \times 10 \text{ cm}^2$
- Cadmium pinhole camera 0.5, 1 and 2 mm in diameter with adjustable height
- Magnification factor $\sim 0.5\text{-}2 \times$
- Position sensitive GEM detector

In current work we are studying in more detailed way the features of our system

- A pinhole camera properties
- Position resolution of the overall system (mainly a pinhole camera diameter and height)
- Pigments recognition enhancement by using Cu anode for the X-ray tube







Main parameters of the GEMROC ASIC

- 32 channels, each split into energy and timing sub-channel
- Self triggering mode
- Data derandomization and zero suppression
- Equipped with testability functions
- Deliver energy, position and time informations

GEM detector parameters

- Position sensitive detector with active area $10 \times 10 \text{ cm}^2$
- 128×128 orthogonal readout strips
- Charge collected on several readout strips
- Position reconstruction
 - ▶ Matching time stamps of signals of X and Y strips
 - ▶ Fine position reconstruction using signals amplitudes (Centre of Gravity method)
- Ar/CO₂ active gas mixture
 - ▶ Low cost
 - ▶ Moderate detection efficiency





Parameters of the system

- Dedicated readout electronics
 - ▶ Simultaneous measurement of time and amplitude
 - ▶ Finding coincidences of X and Y signals (reconstruction of 2D position)
- Energy resolution of about 20 % FWHM for 5.9 keV
- Position resolution of about 100 μm rms. (without pinhole camera)
- High count rate capability up to 5×10^6 cps

Photograph of the whole phantom


Phantom prepared by National Museum in Krakow



Prepared phantom

- A pattern of stripes painted with inorganic pigments on a wooden panel
- The phantom was prepared according to the XV century painting techniques and using paints based on historical pigments by Kremer Pigmente GmbH & Co. KG.
- Different pigments are overlapping each other and are painted on top of three different bottom layers (wood, lead-white and chalk)

Phantom with stripes (zoomed)




- Investigated area of the phantom

Prepared phantom

- A pattern of stripes painted with inorganic pigments on a wooden panel
- The phantom was prepared according to the XV century painting techniques and using paints based on historical pigments by Kremer Pigmente GmbH & Co. KG.
- Different pigments are overlapping each other and are painted on top of three different bottom layers (wood, lead-white and chalk)

Phantom with stripes (zoomed)




- Investigated area of the phantom

Prepared phantom


- A pattern of stripes painted with inorganic pigments on a wooden panel
- The phantom was prepared according to the XV century painting techniques and using paints based on historical pigments by Kremer Pigmente GmbH & Co. KG.
- Different pigments are overlapping each other and are painted on top of three different bottom layers (wood, lead-white and chalk)

Phantom with stripes (zoomed)




- Schematic view of painting layers

Total intensity map





Cumulative energy spectrum





Cr window: 5.0-5.8 keV

- Visible chrome green bottom layer
- Attenuation by Vermilion Hg pigment clearly seen (top-right corner)





Co and Fe window: 6.0-7.3 keV

- Visible cobalt blue and sienna layers





Co window: 6.8-7.0 keV

- Visible cobalt blue





Fe window: 6.3-6.5 keV

- Obtained by subtraction (Co&Fe window – Co window)
- Visible sienna layer





Cu and Zn window: 7.6-9.0 keV

- Visible malachite and zinc white layers





Zn window: 8.55-8.75 keV

- Obtained by subtraction (Cu&Zn window – Cu window)
- Visible zinc white layer



Cu window: 7.9-8.1 keV

- Obtained by subtraction (Cu&Zn window – Zn window)
- Visible malachite layer



Remarks

- The technique is suitable for qualitative pigment determination
- The limitations are mainly due to
 - ▶ Moderate intrinsic energy resolution of the GEM detector
 - ▶ The pigments chemical compositions can be similar (based on the same or slightly different in Z elements)
- X-ray tube with Cu anode provide better visualization of pigments based on lighter elements than Cu (no excitation of the detector inner components)
- X-ray tube with Mo anode is more suitable when pigments of higher-Z elements are dominating (not presented here)
- Future plans: new front-end electronics and Xe-based gas mixture surely will improve the results (to $\sim 15\%$ FWHM of 5.9 keV)


Position resolution of the system – theoretical model

- Pinhole effective aperture: $d_{eff} = d \left[1 - \frac{\Delta x}{s} \left(1 + \frac{1}{M} \right) \right]$

$$\sigma^2 = \frac{d_{eff}^2}{8\ln(2)} \left(1 + \frac{1}{M} \right)^2 + \left(\frac{p}{\sqrt{12M}} \right)^2 + \frac{(x_d \Delta x / s)^2}{8\ln(2)M^2} \left(1 + \frac{1}{M} \right)^2$$

- Detector intrinsic position resolution
- Prallax error depending on the thickness of the drift gap in GEM

- d – diameter of the pinhole
- d_{eff} – effective diameter of the pinhole depending on the photon angle
- M – magnification of the detection system
- p – pitch of the GEM readout structure
- x_d – thickness of the drift gap in GEM
- s – object to detector window distance
- Δx – distance from the centre of the detector to the photon registration position






AGH

Position resolution of the system – theoretical model

- Pinhole effective aperture: $d_{eff} = d \left[1 - \frac{\Delta x}{d} \left(1 + \frac{1}{M} \right) \right]$



$d = 1 \text{ mm}$ and $M = 1$





Δx – distance from the centre of the detector to the photon registration position

Δx



Pinhole diameter $d=2$ mm





Pinhole diameter $d=1$ mm



Pinhole diameter $d=0.5$ mm





Results for different magnification factors $M = \frac{DP}{PS}$



Stainless steel sheet

- Size of plate: 9 × 5 cm
- Pinhole d = 0.5 mm
- Magnification M ~ 1.28
- Milled with patterns with different opening sizes



Conclusions

- The technique is suitable for fast course screening of large area paintings
- Due to limited energy resolution of the GEM detector it is not suitable for detailed elemental analysis
- The system is able to perform qualitative pigment determination
- Position resolution is at the level of 0.5 mm rms.
- Time of the measurement can be optimized at the expense of position resolution and vice versa
- The overall setup can be easily rearranged for traditional radiography imaging (not presented here)

Acknowledgements

- We thank the RD51 collaboration (especially L. Ropelewski) for its support and providing us with the GEM detector.
- Laboratory of Analysis and Non-Destructive Investigation of Heritage Objects, National Museum in Krakow (Ł. Bratasz)

Thank you for your attention!